

#### Overview

High throughput screening (HTS) for neurotoxicity has lagged behind other applications due to the lack of physiologically relevant human neuronal cell models. A co-culture system comprising of human neurons & Schwann cells has remained elusive, even though the advantage of using such a system is well understood. In this study, we have characterized a human iPSCs-derived sensory neurons (hSNs) and primary human Schwann cells (hSCs) co-culture capable of robustly predicting the safety of chemotherapeutic drugs.

### Methods






Plate cells in multiwell plates for high throughput screening



High throughput imaging and analysis



- Neurite length
- Number of branches
- Number of processes Total number of cells

Multielectrode Arrays

- Number of spikes/second
- Number of bursts/second
- Number of spikes/burst
- Interburst interval

Assessing Chemotherapeutics Neurotoxicity Using Human Cells Based High-Throughput Neuronal-Schwann Cells Culture System

Sharma AD<sup>1</sup>, McCoy L<sup>1</sup>, Jacobs E<sup>1</sup>, Curley JL<sup>1</sup>, Moore MJ<sup>1,2</sup> <sup>1</sup>AxoSim, Inc, New Orleans, LA

<sup>2</sup>Biomedical Engineering, Tulane University, New Orleans, LA

# Results – High Throughput Imaging




Fig.1 Images showing monocultures and cocultures exposed to 10µM of Paclitaxel and Oxaliplatin.




Table 1. EC50 and LD50 comparison for mono and co-cultures.

	Paclitaxel		Oxaliplatin	
	EC50-Coculture (uM)	EC50-Monoculture (uM)	EC50-Coculture (uM)	EC50-Monoculture (uM)
Neurite length	0.92	0.23	40.06	23.37
Total processes	5.41	0.69	43.54	20.07
Total branches	0.07	0.09	17.98	11.97
Neuron count	≥100	1.95	42.2	34
Schwann cell count	0.04		19.79	

AxoSim, Inc, New Orleans, LA Correspondence: lowry.curley@axosim.com

# Results – Multielectrode array




Fig.1 Comparison of spontaneous firing rate of cocultures and monocultures. Cocultures showed higher firing rate as compared to monocultures.

# Summary & Conclusion

- Co-cultures and monocultures displayed differences in dose response for both Paclitaxel and Oxaliplatin.
- Neurite outgrowth parameters displayed higher toxicity a compared to cell counts.
- Effective concentration (EC50) for co-cultures was higher than monocultures for both drugs.
- Paclitaxel showed a toxic dose response with Schwann cells as compared to neurons.
- Oxaliplatin did neurotoxicity show not monocultures as compared to co-cultures.
- 30K Peri4U's with 7.5K SCs were active by Day 6.
- Neurons-only condition was slowest to become electrically active. Also, it only happened for high density conditions such as 70K neurons.

### Acknowledgements

JLC and MJM are co-founders of AxoSim Technologies, LLC. MJM is an associate professor at Tulane University. Funding provided by NIH Phase II STTR Grant (R42TR001270).