

Overview

Antibody-drug conjugates (ADCs) are intended for targeted delivery of highly potent payloads to cancer cells and may cause peripheral neuropathy (PN) by a variety of mechanisms including bystander effect (BE).

In this study, we have characterized the co-culture of human iPSC-derived sensory neurons and primary human Schwann cells (HSCs) using two well-characterized chemotherapeutics. The system was then used to compare the effects of non-targeted antibody drug conjugates (ADCs) and small molecules, using monomethyl auristatin E (MMAE) and MMAF. These experiments show the potential of high content screening (HCS) for investigating mechanisms of PN and assessing new therapeutics.



Methods





Human Schwann cells

Plate as mono-cultures and co-cultures in multi-well plates for HCS



High Content Imaging and Analysis



Development of a High-Content Human Co-culture Model to Investigate Chemotherapy-Induced Peripheral Neuropathy

Laurie McCoy¹, Esther Trueblood², Brian Schimpf², Anup D. Sharma¹ ¹AxoSim, Inc, New Orleans, LA ²Seattle Genetics, Seattle, WA.









Mono-cultures

Fig.1 Images showing mono-cultures and co-cultures exposed to 10µM of paclitaxel and oxaliplatin. **Morphological parameters Cell counts**



Fig.2 Mono-cultures and co-cultures exposed to Paclitaxel and Oxaliplatin. Paclitaxel showed a toxic dose response with Schwann cells as compared to neurons. Oxaliplatin did not show neurotoxicity in monocultures as compared to co-cultures.

Table 1. X50 comparison for mono- and co-cultures treated with known PN-causing chemotherapeutics.

| | Pacli | taxel | Oxaliplatin | | | |
|--------------------|--------------------|----------------------|--------------------|----------------------|--|--|
| | Co-culture (uM) | Mono-culture (uM) | Co-culture (uM) | Mono-culture (uM) | | |
| Neurite length | 0.063 | 0.2669 | 21.4 | 30.38 | | |
| Total processes | >100 | >100 | 20.25 | 53.35 | | |
| Total branches | 0.023 | 0.065 | 6.457 | 12.57 | | |
| Neuron count | >100 | >100 | 15.96 | >100 | | |
| Schwann cell count | 0.055 | | 18.72 | | | |

OUTPUTS

Neurite length Number of branches • Number of processes Total number of cells









Fig.3 Mono-cultures and co-cultures exposed to small molecule payloads and non-targeted ADCs with different physiochemical properties. MMAE ADC shows greater potency in the presence of Schwann cells, presumably due to bystander effect.

Table 2. X50 comparison for mono- and co-cultures treated with small molecule payloads and non-targeted ADCs. High bystander effect MMAE ADC (high BE) shows greater potency in the presence of Schwann cells, whereas MMAF ADC (low BE) does not.

| | Small molecule | | | | Antibody Drug Conjugates | | | |
|-------------------|--------------------|----------------------|--------------------|----------------------|--------------------------|-------------------------|-----------------------|-------------------------|
| | Cys-mc-MMAF | | MMAE | | Cys-mc-MMAF | | MMAE | |
| | Co-culture (nM) | Mono-Culture (nM) | Co-culture (nM) | Mono-culture (nM) | Co-culture (µg/mL) | Mono-culture (µg/mL) | Co-culture (µg/mL) | Mono-culture (µg/mL) |
| Neurite length | >100 | >100 | 0.53 | 0.58 | 568 | 656 | 19 | 209 |

Summary & Conclusions

- A sensitive high content assay for evaluating peripheral neuropathy potential was established using co-cultures of human Schwann cells and iPSC-derived sensory neurons. • Differences between responses in co-culture compared to
- mono-culture were seen during characterization with Paclitaxel and Oxaliplatin
 - in co-culture compared to mono-culture.
 - Paclitaxel showed a lower X50 value for neurite length
 - in co-culture vs mono-culture
 - Oxaliplatin showed a lower X50 value for neuron count
 - Paclitaxel and oxaliplatin both showed dose responsive toxicity to Schwann cells in co-culture.
- In both co-culture and mono-culture, MMAE showed lower X50 than MMAF as a small molecule or an ADC.
- The lower X50 of the MMAE ADC in coculture vs mono-culture indicates that bystander effect may contribute to the mechanisms of ADC-induced peripheral neuropathy.

SeattleGenetics®

Results: Small Molecules and ADCs